PRODIGY Pocket[®]

Blood Glucose Monitoring System

Owner's Manual

PRODIGY Pocket® Blood Glucose Meter

Dear Prodigy Pocket® Meter Owner:

Thank you for choosing the Prodigy Pocket® meter as your blood glucose monitoring system. Please read this manual carefully as it contains important information about your new Prodigy® system. A warranty registration card is included with your system. **Please return the completed card to us or complete warranty online at www.prodigymeter.com.**

Prodigy® meters are designed to help you and your healthcare professionals monitor your blood glucose levels. This owner's manual will help you learn how to use your Prodigy® meter effectively. Customer Care is available 24 hours, 7 days a week by calling toll free **1.800.243.2636**. All questions concerning test results should be directed to your healthcare professional.

Prodigy® meters have the latest technology for blood glucose monitoring and are made with quality components. All Prodigy® meters are easy to use, give you fast, accurate test results with a minimal sample of blood, have large, easy-to-read display screens, and are lightweight and portable for your convenience.

Prodigy® meters:

- **Require No Coding**, allowing you to save time and avoid human error due to improper coding.
- Allow you to perform Alternate Site Testing (AST).
- Have memory and data management capabilities. Prodigy's® **free software** gives you and your healthcare professionals powerful graphic tools to manage your diabetes.

Table of Contents

Letter to Meter Owner.....	2-3
Table of Contents	4-5
Important Safety Instructions	6-8
Important Health-Related Information	9-10
Alternate Site Testing (AST).....	11-13
About Prodigy® Test Strips	14-16
Key Functions of the Test Strip	14
Important Test Strip Information.....	15-16
Introduction to Prodigy® Meters	17-18
Intended Use.....	17
Test Principle	17-18
Contents of the Prodigy Pocket® Meter	19
Setting up the Prodigy Pocket® Meter	20-27
Key Functions of the Meter.....	20
Meter Display	21-22
Setup Steps.....	23-27
Control Solution Testing.....	28-35
Performing a Control Solution Test	31-35
Performing a Blood Glucose Test	36-43
Preparing the Lancing Device.....	36-38
Getting a Blood Sample	38-40
Applying Blood to the Test Strip.....	41-43
Using the Meter Memory.....	44-46
Specifications.....	47

Viewing Results on a Personal Computer	48
Caring for your Prodigy® Meter.....	49-52
Cleaning your System	49
Lancing Device	49
Disinfecting your System	49-50
Storing your System	50
Battery	51
Low Battery	51
Battery Replacement.....	52
System Troubleshooting	53-59
Special Messages	53-54
Error Messages	54-55
Troubleshooting Guide.....	56-59
Information About your Prodigy® Meter.....	60-68
Comparing Meter and Laboratory Results.....	60-62
Performance Characteristics	62
Expected Test Results.....	63
Unexpected Test Results.....	64-65
Accuracy of Meter Measurements	66-68
Symbols Information	69
Warranty Information	70
Limited Lifetime Warranty	70

Important Safety Instructions

Read this before using your Prodigy® meter. The following basic safety precautions should always be taken.

- Your Prodigy® meter, test strips, lancets and lancing device are for single patient use only. Do not share them with anyone, including other family members. Do not use on multiple patients.
- Close supervision is necessary when the device is used by, on, or near children, handicapped persons or invalids.
- Use the device only for the intended use described in this manual.
- Do not use test strips and control solutions with your Prodigy® meter that are not supplied by Prodigy®.
- Do not use the device if it is not working properly, or if it has suffered any damage.
- Before using any product to test your blood glucose, read all instructions thoroughly and practice the test. Do all quality control checks as directed and consult with a diabetes healthcare professional.

- Keep the test strip vial away from children; the vial cap and the test strips can be a potential choking hazard.
- Never chew or swallow a test strip. If this occurs, please seek medical assistance immediately.
- All parts of this kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection.
- These devices are intended to be used for patient self-monitoring and should not be used to collect blood from more than one person as this poses a risk of transmitting blood-borne pathogens such as Hepatitis B or HIV.

For further information, please refer to:

“FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Blood-borne Pathogens: Initial Communication” (2010) <http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>

“CDC Clinical Reminder: Use of Fingertick Devices on More than One Person Poses Risk for Transmitting Blood-borne Pathogens” (2010)
<http://cdc.gov/injectionsafety/Fingertick-DevicesBGM.html>

DO NOT CHANGE YOUR TREATMENT BASED ON A SINGLE RESULT THAT DOES NOT MATCH HOW YOU FEEL OR IF YOU BELIEVE THAT YOUR TEST RESULT COULD BE INCORRECT.

Important Health-Related Information

Severe dehydration and excessive water loss may cause false low results. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.

Under normal blood concentrations and conditions, the following should not significantly affect results:

- Elevated blood triglycerides
- Reducing substances such as uric acid and ascorbic acid
- Acetaminophen
- Dopa
- Methyl dopa
- L-dopa
- Tolbutamide

If you are experiencing symptoms that are not consistent with your blood glucose test results and you have followed all instructions described in this owner's manual, contact your healthcare professional.

Test results below 60 mg/dL (3.3 mmol/L) indicate low blood glucose (hypoglycemia). Test results greater than 240 mg/dL (13.3 mmol/L) indicate high blood glucose (hyperglycemia). If your results are below 60 mg/dL or above 240 mg/dL, repeat the test, and if the results are still below 60 mg/dL (3.3 mmol/L) or above 240 mg/dL (13.3 mmol/L), consult your healthcare professional immediately.

Inaccurate results may occur in severely hypotensive individuals or patients in shock. Inaccurate results may occur for individuals experiencing a hyperglycemic-hyperosmolar state. Please refer to the test strip package insert for additional important information.

Alternate Site Testing (AST)

There are important limitations to Alternate Site Testing (AST).

Please consult your healthcare professional before performing AST.

What is AST?

Alternate Site Testing (AST) means you can use parts of the body other than your fingertips to check your blood glucose levels. Prodigy® meters allow you to test on your palm, forearm, upper arm, calf or thigh. See Figure 1.


Figure 1

What is the Advantage?

Fingertips feel pain more readily because they are full of nerve endings (receptors). At other body sites, nerve endings are not so numerous and you will not feel as much pain as you will experience at the fingertip.

When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels. Capillary blood at the fingertip reflects these changes faster than capillary blood at other sites. Therefore, if you are testing your blood glucose level during or immediately after a meal, physical exercise or stressful event, take the blood sample from your fingertip only.

Use AST only:

- ❖ Two hours or more after your last meal.
- ❖ Two hours or more after taking insulin.
- ❖ Two hours or more after exercise.
- ❖ During a relaxed and calm state.

Do not use AST if:

- ❖ You have reason to believe you have hypoglycemia or hyperglycemia.
- ❖ Your routine glucose results are often fluctuating.
- ❖ You are pregnant.

**To increase accuracy when using AST, rub the puncture site before extracting blood. An alternate site device is required to perform AST. If you do not have an alternate site device but would like to use one, please contact Customer Care at 1.800.243.2636.*

About Prodigy® Test Strips

Key Functions of the Test Strip

All Prodigy® meters measure the amount of sugar in the blood, commonly referred to as blood glucose. Blood is applied to the opening of the absorbent channel on the end of the test strip and is automatically drawn into the test strip.

The test strip consists of the following parts:

Absorbent Channel

Apply a drop of blood and it will be drawn in automatically.

Confirmation Window

Shows whether enough blood has been drawn into the test strip's absorbent channel.

Test Strip Handle

Hold this part to insert the test strip into the test strip port on the meter.

Contact Bars

Insert this end of the test strip into the meter. Push it in firmly until it will not go any further.

Figure 2

Please refer to the *“Performing a Blood Glucose Test”* section for complete instructions.

Important Test Strip Information

- Store test strip packages in a cool, dry place between 39°F–104°F (4°C–40°C). Keep away from direct sunlight and heat. Do not refrigerate. Do not store in your car.
- Store your test strips in their original vial only. Do not transfer them to a new vial or any other container.
- With clean, dry hands you may touch the test strip anywhere on its surface when removing it from the vial or inserting it into the meter.
- Do not clean the testing site with alcohol. Instead, wash the test site thoroughly with soap and water only.
- Immediately use a test strip after removing it from the vial, replace the vial cap and close it tightly to prevent exposure to moisture.
- Only apply a blood sample or a control

solution sample to the test strip's absorbent channel. Applying other substances to the test strip's absorbent channel will cause inaccurate results.

- Record the discard date on the vial label when you first open it. **Discard remaining test strips 90 days after the first opening date.**
- Do not use test strips beyond the expiration date printed on the package.
- After removing the test strip from the vial, insert it directly into the meter. Do not lay the test strip on any other surface.

Warning: Keep the test strip vial away from children; the vial cap and the test strips can be a potential choking hazard. Never chew or swallow a test strip. If this occurs, please seek medical assistance immediately.

Introduction to Prodigy[®] Meters

Intended Use

Prodigy[®] blood glucose monitoring systems are intended for use outside the body (in vitro diagnostic use only) and should only be used for blood glucose (blood sugar) testing with fresh capillary whole blood samples.

These systems are intended for use in the home and in clinical settings and should not be used for the diagnosis of diabetes or for the testing of newborns. Consult your healthcare provider if unusual readings occur.

Test Principle

Blood glucose is measured by an electrical current that is produced when a blood sample mixes with the reagent (special chemicals) of the test strip. The electrical current changes with the amount of glucose in the blood sample. Prodigy[®] meters measure the strength of the electrical current, calculate your blood glucose level and then display your result in either mg/dL or mmol/L.

PRODIGY Pocket® Blood Glucose Meter

Prodigy® meters, test strips and control solutions have been designed, tested and proven to work together as a system to produce accurate blood glucose test results.

Important: Use only Prodigy® control solutions and test strips with your Prodigy® meter. Using other test strips and control solutions with this meter can produce inaccurate results.

Contents of the Prodigy Pocket Meter®

The Prodigy Pocket® meter is available as a meter only or as a meter kit. Please check the “REF” number marked on the outside of the box to see if you have purchased a “Meter” or a “Meter Kit.” Please review the contents of your purchase to confirm that all the components are included as listed on the side of your meter box.

Setting up the Prodigy Pocket® Meter

Key Functions of the Meter

Test Strip Slot

Insert the test strip here. The meter will turn on automatically.

LCD Display

Guides you through the test using symbols and simple messages.

Main Button

Turns the meter on or performs other functions described in this manual.

Set Button

Located on the back of the meter, inside the battery compartment; used to setup the meter.

Data Port

Port for USB cable connection to your computer.

Figure 3

Figure 4

Figure 5

Meter Display

Control Solution Test Symbol

Shows that you are in control solution mode. Your test result will not be stored in meter memory.

Test Result Area

Displays test results.

Time Date

Blood Drop Symbol

Flashes when sample should be applied.

Figure 6

Unit of Measure

Indicates what unit of measure the test result is displayed in.

Meter Display

Temperature Symbol

Appears when temperature is outside of operating range.

Low Battery Symbol

Appears when the battery power is low.

Figure 7

Memory Symbol

Appears when you review the memory.

Test Strip Symbol

Appears when the meter is in testing mode.

Day Average

Indicates that the displayed test result is an average.

Setup Steps

The Prodigy Pocket® meter has several different settings that may need to be adjusted before using the meter. If you replace the battery, you should check and update the time and date.

STEP 1: Insert the Battery

Open the battery cover located on the backside of the meter. Insert one (1) 3 Volt CR2032 lithium battery. The meter should beep, turn on and enter Setting Mode. See Figure 8.


Figure 8

STEP 2: Enter Setting Mode

If your meter is off, press and hold the "SET" button located in the battery compartment. The meter turns on and is now in setting mode.

STEP 3: Set the Year

The year appears with the number flashing. Press and release the "M" button to advance the year.

With the correct year on the display, (See Figure 9) press the **"SET"** button and the date will appear on the display with the month segment flashing. See Figure 10.


Figure 9

STEP 4: Set the Month

Press and release the **"M"** button to advance the month. With the correct month on the display, press the **"SET"** button and the date segment will start flashing. See Figure 11.


Figure 10

STEP 5: Set the Date

Press and release the **"M"** button to advance the date. With the correct date on the display, press the **"SET"** button and the time will appear on the display with the hour segment flashing. See Figure 12.

Figure 11

STEP 6: Set the Hour

Press and release the "M" button to advance the hour. With the correct hour on the display, press the "SET" button and the minutes segment will start flashing. See Figure 13.


Figure 12

STEP 7: Set the Minutes

Press and release the "M" button to advance the minutes. With the correct minute on the display, press the "SET" button and the current unit of measurement will start flashing.


Figure 13

Important: Day averages are calculated from results obtained during the 7, 14 and 28 days preceding the current date and time settings. When the date and time are changed, the 7, 14 and 28-day averages may change.

STEP 8: Unit of Measure

Press and hold the **"M"** button until the unit of measurement you are choosing appears on the display.

See Figures 14 and 15.

Your meter can display test results in milligrams per deciliter (mg/dL) or millimoles per liter (mmol/L).

- The mg/dL unit is standard in the United States.
- The mmol/L unit is standard in Canada.


Figure 14


Figure 15

Press the **"SET"** button and **"dEL"** will be displayed. See Figure 16

STEP 9: Delete Memory

When the **"dEL"** symbol and the flashing **"M"** symbol appear on the display, you can choose to clear the memory. If you do not want to clear the memory, press the **"SET"** button again to skip this step. If you want to clear ALL memory, press and hold the **"M"** button for four (4) seconds.

The " - - - " image will appear on the LCD screen to indicate that all memory has been deleted.

Setup is now completed.

"OFF" is displayed and the meter turns off.

See Figure 17

Figure 16

Figure 17

Control Solution Testing

If your Prodigy® control solution did not come with a Prodigy® Control Solution Insert, the following information shall take the place of the insert.

Product Numbers:

990310 (Low)	53310 (Low)	53350 (High)
--------------	-------------	--------------

Intended Use:

Prodigy® control solution is a red liquid that contains glucose that will react with test strips and produce a test result. You can use either high or low Prodigy® control solution with your Prodigy® system.

Summary:

Prodigy® control solutions contain certain amounts of glucose that can be tested with all Prodigy® brand glucose meters. Test results should fall within the control ranges printed on the back of your Prodigy® test strip vial.

Caution: The control solution ranges are located on the back of the Prodigy® test strip vial. They are not recommended target ranges for your blood glucose.

Important Information:

- Use only Prodigy® control solutions with your Prodigy® meter.
- Check the expiration date on the bottle. Do not use if expired.
- Use within a period of 90 days from the date that you first open it. Record the discard date on the control solution bottle when you first open it to serve as a reminder to discard after 90 days.
- For in vitro diagnostic use.

Why Perform a Control Solution Test:

- To ensure that your meter and test strips are working properly together.
- To allow you to practice testing without using your own blood.

When to use:

- Once a week (to ensure that you continue to have accurate results).

PRODIGY Pocket® Blood Glucose Meter

- When you begin using a new vial of test strips.
- When test strips have been exposed to extreme environmental conditions.
- If you drop the meter.
- If you change the battery in the meter.

Storage and Handling:

- Store the control solution tightly closed at temperatures below 86°F (30°C).
- Do not refrigerate. Do not store in your car.

Composition:

	Glucose (% w/v)	Nonreactive Ingredients (% w/v)
Low Solution	0.05	99.95
High Solution	0.19	99.81

Performing a Control Solution Test

Start with the meter off.

STEP 1: Wash your Hands

Wash your hands with mild soap and water before performing any test. Be sure to dry them thoroughly.

STEP 2: Insert Test Strip

Insert a test strip with the contact bar end entering into the test strip slot first. Push the test strip as far as it will go without bending it. The meter turns on automatically.

STEP 3: Mark as a Control Solution Test

After the "💧" symbol appears on the display, press the "M" button and a "🧴" symbol appears on the display. With the "🧴" symbol on the display, the meter will not store your test result in the memory. If you decide not to perform a control solution test, press the "M" button again and the "🧴" symbol will disappear.

Important: Be sure that you are in control solution mode so that the test result will not be stored in the meter memory.

STEP 4: Apply Control Solution

1. Check the expiration and discard dates on your control solution and test strip vials.
2. Shake the control solution bottle well, then remove cap.
3. Squeeze the bottle and discard the first drop then wipe the dispenser tip with clean tissue paper or cotton.
4. Squeeze the bottle again to get a second drop onto a clean, non-absorbent surface or on your clean fingertip. See Figure 18.
5. Bring the tip of the test strip to touch the drop of solution until the meter beeps.


Figure 18

Caution: To avoid contaminating the control solution with the content of the test strip, **DO NOT DIRECTLY APPLY CONTROL SOLUTION ONTO THE TEST STRIP.** See Figures 19 and 20.

Figure 19


Figure 20

STEP 5 : Check if the Test Result is in Range.

After the meter counts down from 6 to 1, the control solution test result appears on the display.

Compare the test result with the range printed on the test strip vial. Each vial of Prodigy® No Coding Test Strips may have a different control solution range. The result should fall within the printed range on the test strip vial. See Figure 21.


Figure 21

Questionable Results:

If the test result falls outside the range printed on the test strip vial, check the "Trouble-Shooting Guide" located in the "System Troubleshooting" section of this manual and repeat the test.

Results that fall outside of the expected control ranges may indicate:

- An error in performing the test
- Control solution bottle not shaken well
- Expired or contaminated control solution
- Control solution that is too warm or too cold
- Test strip deterioration
- Meter malfunction
- Test strips expired

Caution: Do not use test strips or control solution that have exceeded the discard date, are expired or have been damaged. Your results may be inaccurate.

Contact Customer Care

If your test strips have been damaged or if you continue to get out-of-range results, it means that the system or the control solution may not be working properly. **DO NOT** use the system to test your blood glucose level. If you are unable to resolve the problem, contact Customer Care at **1.800.243.2636**.

Performing a Blood Glucose Test

Preparing the Lancing Device

STEP 1: Remove the cap of the lancing device by twisting it off. See Figure 22.


Figure 22

STEP 2: Insert a sterile lancet into the lancet holder of the lancing device and push down firmly until the lancet is fully seated. Do not twist the lancet. See Figure 23.


Figure 23

STEP 3: Remove the protective cap from the lancet by twisting it and then save it for later use. See Figure 24.


Figure 24

STEP 4: Replace the cap onto the lancing device. Screw the cap until it is snug but not too tight.


Figure 25

STEP 5: Set the lancing level. The adjustable tip offers five (5) levels of skin penetration. To select the desired depth, twist the adjustable tip in either direction

until the number lines up with the arrow. Choose from the following depths: 1-2 (least penetration) for soft, 3 (average penetration) for average, 4-5 (deepest penetration) for thick or calloused skin. See Figure 25.

Warning: To reduce the risk of infection:

- Never share a lancet or lancing device.
- Always use a new, sterile lancet.
- Lancets are for single use only, do not reuse.
- Avoid getting hand lotion, oils, dirt or debris in or on the lancets and the lancing device.
- These devices are intended to be used for patient self-monitoring and should not be used to collect blood from more than one person as this poses a risk of transmitting blood-borne pathogens such as Hepatitis B or HIV.

STEP 6: Cock the lancing device. Slide the cocking control back until it clicks. If it does not click, the lancing device may have been cocked when the lancet was inserted. See Figure 26.

Figure 26

*The lancing device is prepared and ready to lance your finger for a blood sample.

Getting a Blood Sample

STEP 1: Wash Your Hands and Puncture Site

Use warm, soapy water. Rinse and dry your hands thoroughly.

STEP 2: Select and Lance a Puncture Site

- **Fingertip**

Hold the lancing device firmly against the side of your finger. Press the release button. You will hear a click, indicating that the puncture is complete. See Figure 27.


Figure 27

- **Sites other Than Your Fingertip**

Please refer to the “About Alternate Site Testing (AST)” Section. Please consult your healthcare professional before obtaining blood from sites other than your fingertip.

STEP 3: Obtain a Blood Sample

Do not smear the blood sample. See Figure 28. To obtain the most accurate results, wipe off the first drop of blood and gently milk the finger to obtain another drop of blood. See Figure 29.

Figure 28


Figure 29

STEP 4: Remove the Lancet

Remove the cap of the lancing device by twisting it off. Place the protective cap back on the exposed tip of the lancet by pushing the lancet into it. Eject the lancet by sliding the eject button forward. See Figure 30.


Figure 30

*Always use caution when removing the lancet. Discard the lancet according to your local safety regulations.

Warning: The first drop of blood usually contains tissue fluid and serum, which may affect the test result. It should be discarded.

Applying Blood to the Strip

STEP 1: Insert the Test Strip

Insert a test strip with the contact bar end entering into the test slot first. Push the test strip as far as it will go without bending it. The meter turns on automatically. See Figure 31.

Figure 31

STEP 2: Apply Blood Sample

When the meter shows the “”, apply blood to the opening of the absorbent channel of the test strip where it meets the narrow channel. Blood will be drawn into the test strip. See Figure 32.

The test strip confirmation window should be completely filled before the meter begins to count down. The meter will beep when the confirmation window is full. See Figure 33. If the confirmation window is not full, “**L-b**” is displayed.

Figure 32

Figure 33

You will have to use a new test strip.

If you have trouble filling the test strip, please contact Customer Care at **1.800.243.2636** for assistance.

Caution:

DO NOT:

- Smear or scrape the blood onto the test strip.
- Apply blood to the test strip when the test strip is out of the meter.
- Put blood or foreign objects into the test strip port.

STEP 3: Read Your Result

After the meter counts down from 6 to 1, your blood glucose test result appears along with the unit of measure, date and time. This blood glucose result is automatically stored in the meter memory. Turn the meter off by removing the test strip. Discard the used test strip carefully to avoid contamination. See Figure 34.

Figure 34

Important: If you do not apply a blood sample within three (3) minutes, the meter will automatically turn off. You must remove the test strip and re-insert it again to turn on the meter and restart the test procedure.

Caution: If you cannot test due to a problem with your testing supplies, contact Customer Care at **1.800.243.2636**. Failure to test could delay treatment and/or lead to serious medical conditions.

Using the Meter Memory

The Prodigy Pocket® stores the 120 most recent blood glucose test results with date and time in its memory. It also provides you with 7, 14 and 28-day averages of your blood glucose test results. You can review the individual or average test results by entering the memory mode.

STEP 1: Enter the Memory Mode

While the meter is turned off, press the "M" button to turn the meter on. Press the "M" button again to enter memory mode. The 7-day average

Figure 35

will appear, indicating that you are in the memory mode. If you continue to press the "M" button, the 14 and 28-day averages will appear in order. You can then review the last 120 individual test results in the memory. When using the meter for the first time, " - - - " appears, showing that there are no test results in the memory. See Figure 35.

STEP 2: Recalling Average Test Results

The 7-day average is calculated from the blood glucose results obtained during the last 7 days. It also indicates how many blood glucose tests have been performed within this period, e.g., 21 (21 tests in the last 7 days). See Figure 36.

Figure 36

The 14-day average is calculated from the blood glucose results obtained during the last 14 days. See Figure 37.

Figure 37

It, too, indicates how many blood glucose tests have been performed, e.g., 41 (41 tests in the last 14 days). The 28-day average shows the same information.

STEP 3: Recalling Individual Test Results

After the 28-day average, the most recent test result with date and time will be shown. Press the "M" button once and the next most recent test result will appear.

Each time you press and release the “M” button, the meter will recall up to your last 120 test results in order. See Figure 38.

Figure 38

When the memory is full, the oldest result is deleted as the newest is added. After reaching the last set of results, press the “M” button and the meter turns off. See Figure 39.

Figure 39

STEP 4: Exit the Memory Mode

Press and hold the “M” button for three (3) seconds to turn off the meter at any point.

Important: If you do not press any button for one (1) minute, the meter will show “OFF” and turn off automatically.

PRODIGY Pocket® Blood Glucose Meter

Specifications

Dimensions: 3.54 in. (L) x 1.38 in. (W) x 0.83 in. (H)
90 mm (L) x 35 mm (W) x 21 mm (H)

Weight: 1.34 oz with battery, 38 g

Power Source: One (1) CR2032 lithium battery.

Battery Life: Over 1,000 tests

Display: LCD

Memory: 120 test results with date and time

External Output: USB Data Port

Auto electrode inserting detection

Auto sample loading detection

Auto reaction time count-down

Temperature warning

Operating Condition: 50°F - 104°F (10°C - 40°C),
below 80% R.H. (non-condensing)

Storage/Transportation Condition:
39.2°F - 104°F (4°C - 40°C), below 80% R.H.

Measurement Units: mg/dL or mmol/L

Measurement Range: 20 - 600 mg/dL (1.1 - 33.3 mmol/L)

The specifications may be changed without prior notice.

Viewing Results on a Personal Computer

Test results in your meter's memory can be transmitted to your personal computer. Please visit the Prodigy® website, www.prodigymeter.com for the latest instructions.

Caring for your Prodigy® Meter

Cleaning your System

Your Prodigy® meter should be cleaned whenever it is visibly dirty by wiping the outside of the meter using a cloth dampened with either mild detergent mixed with water or 70% isopropyl (rubbing) alcohol. Do not use bleach or other harsh abrasives to clean your meter.

Do not get any liquids, dirt, dust, blood or control solution inside the meter through the test strip port or the data port. Never spray cleaning solution on the meter or immerse it in any liquid.

Lancing Device

Wipe your Prodigy® lancing device with a soft cloth dampened with water and mild detergent or 70% isopropyl (rubbing) alcohol. Do not immerse the lancing device in any liquid. Do not use bleach.

Disinfecting your System

If your Prodigy® meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing

device should be decontaminated prior to use by the second person.

Disinfect your Prodigy® lancing device once a week. After cleaning the device, unscrew the cap and place it in 70% isopropyl (rubbing) alcohol for several minutes. Make sure the cap is completely dry before reattaching it to the device.

Please contact Customer Care if you have any questions at **1.800.243.2636**.

Storing Your System

Prodigy® meters do not require special maintenance.

- Store the meter, test strips and control solution in the carrying case after each use in a cool, dry place.
- Do not refrigerate.
- Do not store in the kitchen or bathroom.

*Your meter is a precision instrument. Please handle it with care.

Battery

The Prodigy Pocket® meter comes with one (1) CR2032 lithium battery. Batteries are manufactured in a charged state and not designed for recharging. Recharging the battery can cause battery leakage, or in some cases, high pressure rupture.

Warning: Batteries can explode or leak and cause burns if installed backwards, disassembled, charged or exposed to water, fire or high temperatures.

Low Battery

The meter will alert you when the power is getting low by displaying two (2) different messages:

1. When the "" symbol appears alone on the display, the meter is functional and the results remain accurate, but you should change the battery as soon as possible.
2. When the "" symbol appears together with the "**E-b**" symbol on the display, the battery does not have enough power for a test. You must change the battery before using the meter.

Battery Replacement

When replacing the battery, make sure the meter is turned off.

STEP 1: With the meter off, press the buckle on the battery cover and lift up to open the cover.

STEP 2: Remove the old battery and replace with new one. Be sure to align the battery properly.

STEP 3: Close the battery cover. If the meter does not power on after you have replaced the battery, check that the battery is correctly installed.

STEP 4: Turn the meter on by pressing the “M” button to check the time and date. Replacing the battery does not affect the meter’s memory (previous test results stored in the memory). However, the date, time and unit of measurement settings may need to be updated.

Caution: As with all small objects, the battery should be kept away from small children. If the battery is swallowed, seek medical assistance immediately.

System Troubleshooting

Special Messages

Special messages and error messages help to identify certain problems but do not appear in all cases when a problem has occurred. Improper use may cause an inaccurate result without producing an error message.

In the event of a problem, refer to the information under “Action” in the “Error Messages” section. If you continue to have a problem, please refer to the “Troubleshooting Guide” section. If you follow the actions recommended but the problem is not resolved, please contact Customer Care at **1.800.243.2636** for assistance.

Message	What it means
	“Lo” appears when your result is below the measurement limit, which is less than 20 mg/dL (1.1 mmol/L). “Lo” indicates hypoglycemia (low blood glucose). You should immediately consult your healthcare professional.

Message	What it means
	<p>"HI" appears when your result is above the measurement limit, which is higher than 600 mg/dL (33.3 mmol/L). You should immediately consult your healthcare professional.</p>
Error Message	Description
	<p>What it means: Appears when the battery cannot provide enough power for a test.</p> <p>Action: Replace the battery immediately.</p>
	<p>What it means: Appears when inserting a used test strip.</p> <p>Action: Test with a new test strip. If the problem persists, please contact Customer Care at 1.800.243.2636.</p>

Error Message	Description
E-t	<p>What it means: Appears when the temperature is out of the system operating range (50°F-104° or 10°C-40°C).</p> <p>Action: Repeat the test after the meter and test strip are within the operating temperature range.</p>
L-b	<p>What it means: Insufficient blood sample.</p> <p>Action: Insert a new test strip.</p>

Only the most common error messages are listed. If your meter displays an error message that is not listed, please contact Prodigy Customer Care at 1.800.243.2636.

Troubleshooting Guide

- ❖ **The meter does not display a message after inserting a test strip**

Probable Cause	Actions
Battery exhausted.	Replace the battery.
Battery incorrectly installed or absent.	Check that the battery is correctly installed.
Test strip inserted upside down or incompletely.	Insert the test strip correctly with the bar end entering into the test strip port first.
Defective meter.	Please contact Customer Care at 1.800.243.2636 for assistance.

- ❖ **The test does not start after applying the sample.**

Probable Cause	Actions
Insufficient blood sample	Repeat the test using a new test strip with a larger blood sample.

- ❖ The test does not start after applying the sample.

Probable Cause	Actions
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic shut-off. (Two (2) minutes after user action).	Repeat the test with a new test strip. Apply sample only when the "💧" symbol appears on the display.
Defective meter.	Please contact Customer Care at 1.800.243.2636 for assistance.

❖ **The control solution test is out of range.**

Probable Cause	Actions
Error in performing the test.	Read the instructions thoroughly and repeat the test again.
Control solution bottle not shaken well.	Shake the control solution bottle vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiration date and the discard date of the control solution.
Control solution that is too warm or too cold.	Control solution should come to room temperature (less than 86°F/30°C) before testing.
Test strip deterioration.	Please repeat the test with a new test strip.

- ❖ The control solution test is out of range.

Probable Cause	Actions
Meter malfunction.	Please contact Customer Care at 1.800.243.2636 for assistance.

Information About your Prodigy® Meter

Comparing Meter and Laboratory Results

The test results you obtain from your meter may differ somewhat from your laboratory results due to normal variation. Meter results can be affected by factors and conditions that do not affect laboratory results in the same way (see test strip package insert for typical accuracy and precision data as well as important information on limitations). To make an accurate comparison between meter and laboratory results, follow the guidelines below.

Before you go to the laboratory:

- Perform a control solution test to make sure that the meter is working properly.
- It is strongly recommended to fast for at least eight (8) hours before doing comparison tests.
- Take your meter with you to the laboratory.

While at the laboratory:

- Make sure that the samples for both tests (the meter test and the laboratory test) are

taken and tested within 15 minutes of each other.

- Wash your hands before obtaining a blood sample.
- Never use your meter with blood that has been collected in a gray top test tube.
- Use fresh capillary blood only.

You may still have a variation from the result because blood glucose levels can change significantly over short periods, especially if you have recently eaten, exercised, taken medication or experienced stress.¹ In addition, if you have eaten recently, the blood glucose level from a finger stick can be up to 70 mg/dL (3.9 mmol/L) higher than blood drawn from a vein (venous sample) used for a lab test.²

Therefore, it is best to fast for eight (8) hours before doing comparison tests. Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of bodily fluids (severe dehydration) may also cause a meter result to be different from a laboratory result.

References: ¹Surwit, R.S., and Feinglos, M.N.: Diabetes Forecast (1988), April, 49-51.
²Sacks, D.B: "Carbohydrates." Burtis, C.A, and Ashwood, E.R.(ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.

Performance Characteristics

- **Accuracy:** ± 15 mg/dL when glucose < 75 mg/dL
 $\pm 20\%$ when glucose > 75 mg/dL
- **Precision:** This study shows the CV (correlation variation) is less than 5%.
- **The device has been certified to meet the following standards:**
98/79/EC, IEC 60601-1, IEC 61010-1, IEC 60601-1-2, IEC61326, and ISO 15197

Expected Test Results

Please consult your physician to establish your target ranges. The American Diabetes Association (ADA) suggests the following targets for most nonpregnant adults with diabetes. More or less stringent goals may be appropriate for each individual.

Glycemic control Time of Day	ADA suggested reading
Preprandial plasma glucose (before a meal)	70-130 mg/dl
Postprandial plasma glucose (after a meal)	<180 mg/dl

*Source: American Diabetes Association
(2012. <http://www.diabetes.org/living-with-diabetes/treatment-and-care/blood-glucose-control/checking-your-blood-glucose.html>)*

Your Prodigy® meter is designed to help you and your healthcare professional manage your diabetes. You must always rely on your healthcare professional to interpret your test results and to decide how to treat your diabetes.

Unexpected Test Results

Unexpected test results can occur. When this happens, please refer to these cautions.

Low Blood Glucose Test Results

False low results may occur if you are severely dehydrated. If you think you are severely dehydrated, contact your healthcare professional immediately.

If your test result is lower than 70 mg/dL or “Lo” appears on your meter screen, this could be a result of hypoglycemia. Although this result could be due to a testing error, it is best to treat first and retest later. Contact your healthcare professional immediately. See Figure 40.

Figure 40

High Blood Glucose Test Results

If your test result is higher than 180 mg/dL, this could be a result of hyperglycemia. If you are uncertain about this test result, retest. If the result continues to be higher than 180 mg/dL, contact your healthcare professional immediately.

If "HI" appears on your meter screen, this could be a result of severe hyperglycemia (results higher than 600 mg/dL). Retest your blood glucose. If "HI" appears again, contact your healthcare professional immediately. See Figure 41.

Figure 41

Repeated Unexpected Test Results

If you continue to get unexpected test results, check your blood glucose monitoring system by performing a control solution test. Refer to "Performing a Control Solution Test" section in this manual.

Unusual Red Blood Cell Count

A hematocrit range (percentage of red blood cells in the blood) that is extremely high (above 55%) or low (below 30%) can also cause false readings.

Accuracy of Meter Measurements

FAQ: Can two (2) blood glucose meters produce different readings that are from the same blood sample? Does this make one reading more accurate than the other?

Answer: Yes, different meters can produce different readings using the same blood sample.

No, this does not necessarily mean that one result is more accurate than the other.

The reasons for this are as follows:

1. The FDA recognized standard for glucose meter accuracy is ISO 15197: In Vitro Diagnostic Test Systems Requirements for Blood Glucose Monitoring Systems for Self Testing in Managing Diabetes Mellitus. This standard sets accuracy requirements that manufacturers must meet before being approved by the FDA to market their glucose meters.

2. The accuracy requirements for the blood glucose monitoring system—meters and test strips—as stated in ISO 15197 are:

Ninety-five percent (95%) of the individual glucose results shall fall within ± 15 mg/dL (0.83 mmol/L) of the results of the manufacturer's measurement procedure at glucose concentrations < 75 mg/dL (< 4.2 mmol/L) and within $\pm 20\%$ at glucose concentrations ≥ 75 mg/dL (≥ 4.2 mmol/L).

How does this answer the FAQ

The accuracy of the system is determined by a clinical study using 100 blood samples at interval glucose concentrations ranging from < 50 mg/dL to over 400 mg/dL measured using both the manufacturer's meter and a clinical blood glucose analyzer for comparison (reference measurement).

Example

If the manufacturer's reference measurement is 70 mg/dL this means that the meter measurements are considered accurate if they fall within ± 15 mg/dL of 70 mg/dL (from 55 mg/dL to 85 mg/dL).

If meter A produces a measurement of 60 mg/dL and meter B measures 80 mg/dL, then both are in the acceptable range and meet the accuracy requirement.

Each manufacturer has met the System Accuracy requirement as part of its testing to be cleared for marketing. In the System Accuracy, testing a distribution of blood glucose readings will result at the various glucose concentration levels for one manufacturer's system. This is the same for other manufacturer's systems. Because of the range of acceptable readings for the accuracy requirement, identical readings from two different meters on the same blood sample may not happen.

Symbols Information

Symbol	Referent
	Do not re-use. Single use only.
	Consult Operating Instructions.
	Keep away from sunlight.
	Keep dry.
	Temperature limitation.
	Use by.
	Date of manufacture.
LOT	Batch code.
REF	Catalog number.
SN	Serial number.
CONTROL	Control.

Warranty Information

Limited Lifetime Warranty

Prodigy Diabetes Care, LLC extends a limited lifetime warranty to consumers who purchase a new Prodigy® brand meter. Under this limited lifetime warranty, your new meter is covered from the original date of purchase, as long as it has not been modified, altered, or misused. Under this limited lifetime warranty, Prodigy® will replace the meter free of charge if it is defective in material or workmanship.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE MADE. PRODIGY DIABETES CARE, LLC WILL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING DIRECTLY OR INDIRECTLY FROM THE FAILURE OF THE PRODUCT TO PERFORM IN ACCORDANCE WITH SPECIFICATIONS.

Some states do not allow the exclusion or limitation of other express or implied warranties or incidental or consequential damages, so the above limitations or exclusions may not apply to you.

PRODIGY®

Prodigy Diabetes Care, LLC 1.800.366.5901

PO Box 481928
Charlotte, NC 28269

For **Customer Care**
call: **1.800.243.2636**

www.prodigymeter.com

©2013 Prodigy Diabetes Care LLC. Prodigy® and Prodigy Pocket® are registered trademarks of Prodigy Diabetes Care LLC. All other trademarks are property of their respective owners. All rights reserved. POM654 Rev.2 12/13

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